

Cobas Academy Blood Glucose Monitoring

Follow these steps to complete your Cobas Academy Accu-Chek Inform II **online** training:

- 1) Read the PowerPoint, "Orientation on Blood Glucose Monitoring" first, and then follow these remaining steps.
- 2) You will notice when you click on the "Cobas login link", it states "Session Failure". Just continue by clicking on the blue arrow in the upper right hand corner.
- 3) **FIRST LOGIN**
 - a. All instructions for **FIRST LOGIN** are stated on the page. Please follow carefully. Use your Agency/Locum employee number assigned to you. For any issues contact the Point of Care Coordinator, Joseph Stone at 709-454-0510, or email: joseph.stone@lghealth.ca
- 4) **COMPLETE THE CONTENT OF THE COURSE**

Note: If you are having difficulty with this part, log out, and log back in using your employee number **and** password. Sometimes, the training module doesn't show from the initial login.

 - a. Under "My To-Do List", click on the little blue triangle to start the Accu-Chek Inform II training course;
 - b. Click **Launch**;
 - c. Follow the instructions.
 - d. Complete **ALL** modules, and answer **ALL** questions within each module (1, 2, and 3). You can complete the online exam only after you have completed all three modules.
- 5) **EXAM**
 - a. Close out the window with the training course, and highlight **Certification Exam** on the original page;
 - b. Click the blue triangle to start;
 - c. Click **Start**, and then **Next**.
 - d. Read the instructions carefully;
 - e. Complete the exam once it is started, or you risk losing your work;
 - f. Passing requires 18/20 or 90% correct. You will have 3 attempts to pass the exam. If you are not successful after 3 attempts, contact Joseph Stone. Good luck!
- 6) **PRACTICAL COMPONENT**
 - a. When you arrive at Labrador-Grenfell Health, you will complete the practical component of the training. Your observer can be any RN. Follow these instructions:
 - i. Log onto the glucometer and do both a high and low **Quality Control Test** (in any order);
 - ii. Select **Patient Test** and then **Observer Login**;
 - iii. Give the meter to the observer to login, then take it back and perform a patient test;
 - iv. Hand the meter back to the observer so they can answer the two questions (verifying that the test was performed correctly);
 - v. Press the power button to log off;
 - vi. Dock the meter.

Note: To be certified, both the online and practical training must be completed.